

Flake Ice Maker

Applicable area:

Schools, hospitals, research institutes, enterprises and other types of bio-chemical laboratories, hotels, etc.



Model	PLC500
Ice Making Capacity	44lbs/24hrs
Ice Storage Capacity	22 lbs
Consumption	280W
Way of condensation	Air cooling
Way of Influent	Automatic water tap
Tank shell	Stainless steel
Ice shape	Small particals of irregular snow ice
Compressor / Refrigerant	R134a, CFC -free
External Size(W*D*H)	330*480*605mm
Net Weight	70.5 lbs



Model	PLC502
Ice making capacity(kg/24h)	220.5 lbs / 24 hrs
Ice stroage capacity(kg)	88 lbs
Input current(w)	520
Way of condensation	Air cooling
Way of influent	Automatic water tap
Tank shell	Stainless steel
Ice shape	Small particals of irregular snow ice
Compressor / refrigerant	Imports of CFC-free/134a
Overall dimension(L×W×H)	500×611×940mm
Net weight(kg)	133 lbs.
Applicable area	Schools, hospitals, research institutes, enterprises and other types of bio-chemical laboratories, hotels, etc.

Technical Parameter:

Model	PLC501
Ice making capacity(kg/24h)	110 lbs/24 hrs
Ice stroage capacity(kg)	33 lbs
Input current(w)	380
Way of condensation	Air cooling
Way of influent	Automatic water tap
Tank shell	Stainless steel
Ice shape	Small particals of irregular snow ice
Compressor / refrigerant	Imports of CFC-free/134a
Overall dimension(L×W×H)	400×543×720mm
Net weight(kg)	88 lbs
Applicable area	Schools, hospitals, research institutes, enterprises and other types of bio-chemical laboratories, hotels, etc.

