

Water Bath

Features:

1. Internal material: High quality stainless steel.
2. Use microcomputer intellectual temperature control instrument P.I.D, control temperature stably, high precision, digital display, watch directly and clearly.



Control Panel&Display



Drain Valve

Model	PLT200	PLT201	PLT202	PLT203	PLT204	PLT205	PLT206	
Type	Line	1	1	1	2	2	2	
	Hole	1	2	4	6	4	8	
Temperature Range	RT+5 °C ~100 °C							
Temperature	Sensitivity	±0.1 °C						
	Fluctuation	±0.5 °C						
	Precision	PID microcomputer control, touch type, digital display						
	Sensor	Platinum Resistor PT100						
Ambient Temperature	5~30 °C							
Consumption	300W	500W	1000W	1500W	1000W	1500W	2000W	
Power Supply	110/220V±10%,50/60Hz							
Timing Range	0~9999mins							
Net Weight	9 lbs	12.5 lbs	22.5 lbs	28 lbs	14 lbs	17.5 lbs	21 lbs	
Package Size(mm)	360*280*250	480*290*260	790*290*260	1110*330*290	480*420*250	630*420*250	780*420*250	